


Menstrual Cycle: Basics

A menstrual cycle is measured from the first day of your period to the day before your next period.



The average length of a menstrual cycle is 28 days, but everyone's cycle is different. For example, teenagers might have cycles that last 45 days, whereas people in their 20s to 30s might have cycles that last 21 to 38 days.

NOTE: Cycle length and period length are 2 different things.

The menstrual cycle has 4 phases:

1. Menstrual Phase: This is when women have a period. The uterus lining sheds and flows out of the vagina. The period contains blood, mucus and some cells from the lining of the uterus. The average length of a period is 3 to 7 days.
2. Follicular phase: This phase starts on the first day of the period and lasts for 13 to 14 days. Changing hormone levels cause the lining of the uterus to thicken and follicles to grow on the surface of the ovaries. Usually only one follicle will mature into an egg.

3. Ovulation phase: This is when a mature egg is released from an ovary. This usually happens once a month, about 2 weeks before the next period. This is also the time when the chances of getting pregnant are the highest.

4. Luteal phase: In this phase, the egg travels through the fallopian tubes to the uterus. The uterus lining continues to thicken in preparation for pregnancy. In case one gets pregnant, they won't have a period, otherwise they will have a period and the menstrual cycle will start over again. (14 days)

Body Parameters and how they vary during a cycle 📊

FACTOR	MENSTRUAL PHASE	OVULATORY PHASE	LUTEAL PHASE
Core body temperature	Decreases	Starts to increase	Further increases and remains high till next period
HRV	Increases	Further increases	Decreases
Resting HR	Decreases	Inconclusive	Mild increase
HR	Decreases	Increases	Further increases
RR	Decreases	Further decreases	Mild increase
Skin perfusion	Higher	Lower	Inconclusive


MENSTRUAL PHASE:

Core body temperature, RHR and RR ⬇️


HR ⬇️ compared to other phases

HRV, Skin perfusion ⬆️

OVULATORY PHASE:


Core body temperature starts to 

HRV and HR further 

RR and skin perfusion 

RHR inconclusive





LUTEAL PHASE:

Core body temprature, RHR, HR and RR 

HRV 


App Onboarding:

Users are prompted to integrate and import their cycle data from a third party app (if any), input their primary goal for using this powerplug, period length, cycle regularity, cycle length, period dates, health conditions (if any) and any birth control used in the last 6 months.



Integrations

- We only import period start and end dates from the integrated apps.
- As soon as the user integrates the two apps, the data should sync immediately.
- For the users facing issues logging in their Clue app, flag the cases to Terra on this [thread](#).
- Once the cycle data has been imported into UH from the integrated apps, further changes to the data on the third party app won't change the data in UH. Also, any changes made to the data by the user on UH won't be written back to the integrated app.



Cycle tracking: Tracking for basic cycle health


- Input parameters:
 1. **period_ranges**: A list of period ranges, each containing:
 - Start date of a period
 - End date of the same period
 2. **last_cycle_length**: Integer representing the length of the most recent cycle (for irregular cycles)
 3. **typical_cycle_length**: Integer representing the typical cycle length (for regular cycles)

- How algo works:

Take last cycle_start_date > Add typical_cycle_length (for regular cycles) or last_cycle_length (for irregular cycles) to it > You get predicted_cycle_start_date of the next cycle > Take difference of last cycle_end_date and cycle_start_date to get cycle duration > Add to predicted_cycle_start_date to get predicted_cycle_end_date


NOTE: Current algo assumes future cycles will resemble past cycles and doesn't account for external factors that might affect the cycles (like stress, medication, health conditions, etc)

- At any given time, the user would have to be in one of the 6 states:
 1. Menstrual phase: When the user is in their logged menstrual stage (banner shows period start and calculated end dates)
 2. Follicular phase: The phase between menstrual and ovulation phases (banner shows next predicted period dates)
 3. Ovulation phase: When the user is in their ovulation phase (banner shows next predicted period dates)
 4. Luteal phase: Phase between ovulation and predicted menstrual phase dates (banner shows next predicted period dates)
 5. Predicted menstrual phase: Menstrual dates predicted by the algo (banner shows the current predicted dates)
 6. Unconfirmed phase: If user hasn't logged their periods even after 7 days have passed past the predicted menstrual phase dates (no banner shown)







NOTE: At the moment there's a bug in the algo whereby the user's max cycle length gets capped to 30 days once they start using the PP. This isn't an issue for the regular cyclers but for irregular ones, this can result in app showing unconfirmed phase until they log in their next cycle dates.

- The powerplug remains 'in calibration' mode until the following criterias are met:
 1. At least 3 cycle start dates
 2. At least 50% biomarker data for those cycles
 3. BMI ≤ 28 , age is ≤ 40



- App Overview:





Tracking to get pregnant

- Definitions:
 1. Fertile window: 5 days before ovulation day, ovulation day and one day after ovulation day. So, 7 days in total.
 2. Outside fertile window: Any day that's outside the above defined fertile window duration.
 3. Approaching fertile window: The day after the period ends up until the start of the fertile window.
- How algo works:

Divided into 3 phases:

 1. Phase 1: From start of menses upto <40% of cycle duration.

Uses past data (cycle start date+biomarkers) to predict fertile window
 2. Phase 2: Fertile window (>40% to 70% of cycle duration+5 days)

Biomarker algo kicks in and predicts new fertile window (only if enough biomarker data is present)
 3. Phase 3: From 5 days past 70% of cycle duration up until the next period start date

Uses past data to calculate next fertile window





The predicted ovulation window is determined by anchoring to your period start date and detecting a biosignature pattern from your temperature, resting heart rate (RHR), and heart rate variability (HRV). As the cycle advances, additional data improves the algorithm's accuracy in detecting the biosignature. This can sometimes shift the predicted ovulation window, especially because we're forecasting the window in real-time, during the cycle itself.

NOTE (1): Even 1 day of inaccurate data can affect C&O forecasts. However, we use a 3 day rolling window for data prediction, hence, the user will start seeing major shifts in their fertile window predictions if the data remains inaccurate for 3 days continuously.


NOTE(2): The temperature deviation shown in the powerplug tile is different from the skin temperature deviation captured on the home page of the app. This is because we use a 30 day rolling baseline for the calculation of temp deviation for C&O powerplug.

Ref thread [link](#).

- App Overview:



Women’s Health Dashboard



Cycle Tracking Configs Dashboard:

ADMIN /

Cycle Tracking Configs

Id	Email	User Goal	Period Length	Cycle Length	Cycle Regularity	Third Party Apps	Health Conditions	Birth Control Measures
83170	68xnqwt5j@privaterelay.appleid.com	period_monitoring	5	28	regular		none	none
83169	patriciaquilala@gmail.com	period_monitoring	6	31	regular		none	none
83168	96426n7vpn@privaterelay.appleid.com	period_monitoring	3	26	irregular		none	none
83167	kaleighbrooke7@gmail.com	conception	6	37	irregular		thyroid_disorders	none

Filters
USER EMAIL EQ

Filter Clear Filters

This dashboard gives all the details filled by the user while onboarding into their Cycle & Ovulation Powerplug

Menstrual Cycles Dashboard:

ADMIN /

Menstrual Cycles

New Menstrual Cycle

Id	Email	Period Start Date	Period End Date	Perceived Cycle Length	Actual Cycle Length	Predicted Ovulation Day	Calendar Method Ovulation Day	Predicted Next Period Start	Predicted Next Period End	Source	Created At	Updated At	
11309695	Q2hnyashah@gmail.com	Aug 5, 2024	Aug 7, 2024	28		Aug 19, 2024	Aug 19, 2024	Sep 2, 2024	Sep 4, 2024	health_integration	Mar 20, 2025	Mar 20, 2025	View
11309694	ngnvdcpnhrl@privaterelay.appleid.com	Sep 28, 2024	Oct 3, 2024	30	60	Oct 12, 2024	Oct 13, 2024	Oct 28, 2024	Nov 2, 2024	health_integration	Mar 20, 2025	Mar 20, 2025	View
11309693	ngnvdcpnhrl@privaterelay.appleid.com	Jul 29, 2024	Aug 6, 2024	30	60	Aug 13, 2024	Aug 13, 2024	Aug 28, 2024	Sep 5, 2024	health_integration	Mar 20, 2025	Mar 20, 2025	View
11309692	welfgmsuuren@gmail.com	Mar 14, 2025	Mar 16, 2025	40		Mar 30, 2025	Apr 3, 2025	Apr 13, 2025	Apr 15, 2025	health_integration	Mar 20, 2025	Mar 20, 2025	View

Filters
USER EMAIL EQ

Filter Clear Filters

This dashboard shows all of the cycles and their data tracked by the user along with the predicted cycle dates computed by our algorithm.

NOTE: We use this dashboard only to check cycle data and other onboarding details filled by the user. This is intended to give a complete overview of the data being shown to the user on their app.

NOTE: To troubleshoot any technical issue with the user's powerplug, updating their ring's firmware, resets or making them re-onboard the powerplug is not going to help. We would have to flag it to the tech team.

- User feels their ovulation dates predicted by the app aren't correct:

Troubleshooting: Currently its not possible to input/log/correct ovulation day that is predicted by the system. But this is an enhancement planned in the next version for April.

Ref thread [link](#).

- Users complain that their predicted period days aren't correct: Troubleshooting: They can always "log" their period by tapping the calendar icon on the top right next to the edit cycle details button and we'll revise that to recalculate their phases.

Ref thread [link](#).

- The user feels there are days where their temperature has not been recorded (in cycle tracking): When deviation from the recorded baseline is very close to 0, there are no bars shown on the graph. We're working on improvements to the feature so that this is clear going forward. (Currently may also show as "- -" on a given day, on the app)

Ref thread [link](#).

- Users complain that they are seeing the 'in calibration' banner even after logging in past 3 cycles data and wearing the ring continuously: Troubleshooting: We show the "In calibration" if we can't come up with a fertility high confidence. This high confidence is decided on some factors as described.

Ref thread [link](#).

- User notices that predicted ovulation window changes on a given day: We predict the ovulation window for a user based on 3 biomarkers (temperature deviation, RHR and HRV) - hence, if on any day we find that these biomarkers are more (or less) predictive towards a certain day, the algorithm will change basis the stronger suggestive biomarker pattern. Hence, the predicted ovulation window may dynamically change. Ref thread [link](#).
- Users complain that the temperature being shown in their sleep index section and the one being shown inside the powerplug are different.

- Troubleshooting: For the temperature you see on the Sleep Index details page, we calculate it by excluding nap sleep. This gives you an average skin temperature during your main sleep period, without naps affecting the value.

However, the skin temperature in the Cycle Tracking details page is pulled from your baseline temperature, which includes all temp data within sleep boundaries—naps included. Hence, the difference.

Ref thread [link](#).

WIP

- Expanding period logging limits from the current min 3 and max 10 days capability.
- Handling non bleeding users

Non-bleeding users:

"Amenorrhea," or the absence of menstruation,

can be normal during certain periods like before puberty, during pregnancy or breastfeeding, or after menopause, but can also be a sign of various medical conditions or lifestyle factors

.

Here's a breakdown of the different aspects of non-bleeding menstrual periods:

1. What is Amenorrhea?

- Definition: Amenorrhea refers to the absence of menstruation, which is the regular bleeding that occurs during a woman's reproductive years.
- **Types:**
 - **Primary Amenorrhea:** Occurs when a girl who is 16 years old has never had a menstrual period, or if a girl has not developed secondary sexual characteristics (like breast development and pubic hair) by age 13.
 - **Secondary Amenorrhea:** Occurs when a woman who has previously had regular periods stops having them for three months or more.

2. Normal Situations for Non-Bleeding:

- **Before Puberty:** Before a girl starts menstruating, it's normal to not have periods.
- **During Pregnancy:** A missed period is a common early sign of pregnancy.
- **While Breastfeeding:** Exclusively breastfeeding can suppress menstruation.

- **After Menopause:** Menstruation naturally ceases after menopause, which typically occurs between the ages of 45 and 55.
- **Hormonal Contraception:** Some hormonal birth control methods, like the progestogen-only pill, contraceptive injection, or intrauterine system (IUS), can lead to reduced or no bleeding.

3. Potential Causes of Amenorrhea (Other than Normal Situations):

- **Medical Conditions:**
 - **Polycystic Ovary Syndrome (PCOS):** A hormonal disorder that can cause irregular periods or amenorrhea.
 - **Pituitary or Thyroid Disease:** Problems with these glands can disrupt hormone production and affect menstruation.
 - **Eating Disorders:** Conditions like anorexia nervosa can cause amenorrhea due to low body weight and poor nutrition.
 - **Excessive Exercise:** Overexercising, especially in combination with low body weight, can lead to amenorrhea.
 - **Asherman's Syndrome:** A rare condition where scar tissue forms inside the uterus, potentially blocking the flow of menstrual blood.
 - **Genetic Disorders:** Certain genetic conditions, like Fragile X syndrome, can cause menstrual periods to stop early.
- **Lifestyle Factors:**
 - **Stress and Anxiety:** High levels of stress can disrupt the menstrual cycle.
 - **Weight Changes:** Significant weight loss or gain can affect menstruation.
- **Other:**
 - **Travel:** Changes in routine and jet lag can sometimes disrupt periods.
 - **Illness:** Certain illnesses can also cause irregular periods or amenorrhea.

SPOTTING:

Spotting refers to **light vaginal bleeding that occurs outside of your regular menstrual period**

Spotting can have various causes, including hormonal fluctuations, ovulation, birth control methods, pregnancy (implantation bleeding), or even certain medical conditions.

Spotting is lighter and less prolonged than a regular period, and it usually doesn't require the use of pads or tampons

- Make provision to input LH result
- Increase temp data to upto 2 decimal places.
- Freezing ovulation prediction beyond a certain point of time and not letting it keep shifting.

- And many more...

Specific improvements we're working on to address some of the issues over the next few weeks:

1. **Stabilizing the ovulation window** — Once a high-probability ovulation window is identified, we'll prevent multiple shifts later in the same cycle.
2. **Marking ovulation day** — You'll be able to manually mark your ovulation day to refine predictions and improve future accuracy.
3. **Bug fixes for irregular cycles** — We're fixing issues that impact predictions for users with irregular cycle lengths.